

## Contents

COMPETITION ASSAY .....	2
Membrane preparation .....	2
Incubation and filtration .....	2
Data analysis .....	2
SATURATION ASSAY .....	2
Membrane preparation .....	2
Incubation and filtration .....	3
Data analysis .....	3
LIVE CELL ASSAY .....	3
Incubation and filtration .....	3
Data analysis .....	4
PLATE LAYOUTS .....	4
Competition assay .....	4
Saturation assay .....	4

## COMPETITION ASSAY

### Membrane preparation

Frozen tissue or washed cells are partially thawed and homogenized in 20 volumes of cold lysis buffer (50mM Tris-HCl, 5 mM MgCl<sub>2</sub>, 5 mM EDTA, protease inhibitor cocktail). After a low speed spin (1,000 x g for 3 min) to remove large tissue chunks, the homogenate is centrifuged at 20,000 x g for 10 minutes at 4°C to pellet the membranes. The pellet is resuspended in fresh buffer and centrifuged at the same speed a second time. The pellet is then resuspended into buffer (15 ml) containing 10 % sucrose as a cryoprotectant, divided into 1 ml aliquots and stored at -80 °C. A sample of the homogenate is analyzed for protein content using the Pierce® BCA assay. On the day of the assay the membrane preparation is thawed and the pellet resuspended in final assay binding buffer (50 mM Tris, 5 mM MgCl<sub>2</sub>, 0.1 mM EDTA, pH 7.4).

### Incubation and filtration

The filtration binding assay is carried out in 96-well plates in a final volume of 250 µL per well. To each well is added 150 µL membranes (3 - 20 µg protein for cells or 50 - 120 µg protein for tissue), 50 µL of the competing test compound and 50 µL of radioligand solution in buffer. The plate is incubated at 30 °C for 60 minutes with gentle agitation. The incubation is stopped by vacuum filtration onto 0.3 % PEI presoaked GF/C filters using a 96-well FilterMate™ harvester followed by four washes with ice-cold wash buffer. Filters are then dried for 30 minutes at 50°C. The filter is sealed in polyethylene, scintillation cocktail (Betaplate Scint; PerkinElmer) added and the radioactivity counted in a Wallac® TriLux 1450 MicroBeta counter.

### Data analysis

For each drug concentration, non-specific binding is subtracted from total binding to give specific binding. Data is fitted using the non-linear curve fitting routines in Prism® (Graphpad Software Inc). For competition assays, K<sub>i</sub> values are calculated from IC<sub>50</sub> values using the formula  $K_i = IC_{50} / (1 + ([S]/K_d))$  where [S] is the radiotracer concentration used in the assay and K<sub>d</sub> is the dissociation constant of the radiotracer.

## SATURATION ASSAY

### Membrane preparation

Frozen tissue or washed cells are partially thawed and homogenized in 20 volumes of cold lysis buffer (50mM Tris-HCl, 5 mM MgCl<sub>2</sub>, 5 mM EDTA, protease inhibitor cocktail). After a low speed spin (1,000 x g for 3 min) to remove large tissue chunks, the homogenate is centrifuged at 20,000 x g for 10 minutes at 4°C to pellet the membranes. The pellet is resuspended in fresh

buffer and centrifuged at the same speed a second time. The pellet is then resuspended into buffer (15 ml) containing 10 % sucrose as a cryoprotectant, divided into 1 ml aliquots and stored at -80 °C. A sample of the homogenate is analyzed for protein content using the Pierce® BCA assay. On the day of the assay the membrane preparation is thawed and the pellet resuspended in final assay binding buffer.

### Incubation and filtration

The filtration binding assay is carried out in 96-well plates in a final volume of 250 µL per well. To each well is added 150 µL membranes (3 - 20 µg protein for cells; 50 - 120 µg protein for tissue), 50 µL of the unlabeled compound (non-specifics) or buffer and 50 µL of radioligand solution in binding buffer. The radioligand is added at up to 8 different concentrations (e.g. 0.2 - 20 nM). The plate is incubated at 30 °C for 60 minutes with gentle agitation. The incubation is stopped by vacuum filtration onto 0.3 % PEI presoaked GF/C filters using a 96-well FilterMate™ harvester followed by four washes with ice-cold wash buffer. Filters are then dried for 30 minutes at 50°C. The filter is sealed in polyethylene, scintillation cocktail (Betaplate Scint; PerkinElmer) added and the radioactivity counted in a Wallac® TriLux 1450 MicroBeta counter.

### Data analysis

For each radioligand concentration, non-specific binding is subtracted from total binding to give specific binding. Bound CPM values are converted to fmoles per mg protein. Data is fitted using the saturation analysis non-linear curve fitting routines in Prism® (Graphpad Software Inc). The  $K_d$  (in nM) and  $B_{max}$  (fmol/mg or sites/cell) are derived from the saturation curve.

## LIVE CELL ASSAY

### Incubation and filtration

At room temperature, 150 µL of cells (50,000 - 150,000 cells) is added to each well, followed by 50 µL of "cold" unlabeled antibody or protein solution (for non-specifics) or buffer and 50 µL of [<sup>125</sup>I]-labeled antibody or protein. The [<sup>125</sup>I]-labeled antibody or protein is added at up to 8 different concentrations (e.g. 0.13 - 16 nM). The plate is incubated for 120 min with gentle agitation on an orbital shaker. The incubation is stopped by mild vacuum filtration onto GF/C filtermats (presoaked in buffer with BSA) using a 96-well FilterMate™ harvester. Washing is accomplished using an initial wash with ice-cold assay buffer followed by several washes with ice cold PBS. Filters are then dried under a stream of warm air. Scintillation cocktail (Betaplate Scint; PerkinElmer) is added to the filters and the radioactivity counted in a Wallac® TriLux 1450 MicroBeta counter.

## Data analysis

Bound CPM values are converted to fmoles bound [<sup>125</sup>I]protein or antibody per assay well and plotted against [<sup>125</sup>I]protein concentration (nM). Data is fitted using the saturation analysis non-linear curve fitting routines in Prism® (Graphpad Software Inc) to derive the K<sub>d</sub> (in nM) and receptor or surface antigen density (B<sub>max</sub>, fmol/mg or sites/cell) values.

## PLATE LAYOUTS

### Competition assay

Typical plate layout for 4 competing drugs with K<sub>i</sub> values in the range of 3 - 30 nM. Duplicate replicates.

		Competing Drug Concentration												
		0	0.1 nM	0.3 nM	1 nM	3 nM	10 nM	30 nM	100 nM	300 nM	1,000 nM	3,000 nM	Non-specific	
A	Compound 1													
B	Compound 1													
C	Compound 2													
D	Compound 2													
E	Compound 3													
F	Compound 3													
G	Compound 4													
H	Compound 4													

### Saturation assay

Typical plate layout for 2 radiolabeled compounds, each with K<sub>d</sub> values in the range of 0.5 - 2 nM. Triplicate replicates.

		Compound 1			Compound 2									
		Total	Non-Specific		Total	Non-Specific								
		1	2	3	4	5	6	7	8	9	10	11	12	
A	0.13 nM													
B	0.25 nM													
C	0.5 nM													
D	1 nM													
E	2 nM													
F	4 nM													
G	8 nM													
H	16 nM													