


## KINETIC AFFINITY

Kinetic affinity assays provide a technique for determining the affinity of a ligand for its receptor by measuring the real-time binding association and dissociation rates using Surface Plasmon Resonance (SPR). In this method, one of the binding partners (e.g. the receptor protein) is immobilised onto a biosensor surface. The second partner (e.g. the drug ligand) is then continuously flowed across the biosensor surface, where it binds to the immobilised receptor. Binding is measured as an increase in resonance units (RUs) on the biosensor surface. Measuring the increase in binding over time for a given ligand concentration gives the association rate ( $k_a$  or  $k_{on}$ ). Ceasing to flow drug ligand and changing to buffer alone then allows the ligand to wash off the receptor. Measuring the decrease in bound ligand over time gives the ligand dissociation rate ( $k_d$  or  $k_{off}$ ). The affinity of the ligand for the receptor (the equilibrium dissociation constant,  $K_D$ ) is calculated from the kinetic association and dissociation rates ( $k_d/k_a$ ) for several different ligand concentrations.




Kinetic affinity measurements can be complementary to radioligand equilibrium binding affinity, providing additional binding dynamics. For instance, two ligands may have the same binding affinity, but kinetic measurements can show differences in dissociation rates that may lead to selection of a ligand with a slower off-rate as a more effective receptor blocker.

## EQUILIBRIUM BINDING

Receptor-ligand kinetic affinity determinations are usually successful for receptors with larger ligands that have slower association and dissociation rates, e.g. chemokine receptors and cytokine ligands.

For receptor-ligand interactions with faster association and dissociation rates, kinetic affinity analysis can be inappropriate. For these situations, equilibrium binding affinity is more appropriate.



The example shown is for the equilibrium binding of the endogenous estrogen, 17 $\beta$ -estradiol, to human estrogen receptor  $\alpha$  (ER $\alpha$ ). Recombinant his-tagged human ER $\alpha$  was captured onto a nickel-treated NTA biosensor chip surface, then 17 $\beta$ -estradiol from 0.1 to 30 nM was bound to the ER $\alpha$  for 60 seconds, followed by 180 seconds dissociation (only part is shown). Blank-subtracted average binding responses were fitted to a saturation equilibrium binding model to generate a  $K_D$  of 0.18 nM. (IUPHAR reference  $K_D$  for 17 $\beta$ -estradiol binding to ER $\alpha$  = 0.16 nM).

#### LABEL-FREE AFFINITY MEASUREMENT

Surface plasmon resonance is label-free; it does not require either of the binding partners to be labelled with a radioisotopic or fluorescent label. It therefore provides an alternative to a radioligand binding assay where a [ $^3$ H]-radioligand is unavailable or custom radiolabelling is impractical. It also avoids problems on the occasions when incorporation of a bulky isotope (e.g. [ $^{125}$ I)]) might adversely affecting the ligand structure and hence its affinity for the receptor.