

# Techniques for Studying Receptor – Ligand Interactions

Service	Description	Recommended for	Sample / Materials	Deliverables	Pros	Cons
<b>Radioligand Binding</b>	Radiometric technique used to measure ligand-receptor binding parameters (96-well plates)	<ul style="list-style-type: none"> <li>Determining ligand-receptor affinity and receptor expression levels (saturation assays)</li> <li>Determining affinity of unlabelled test compounds competing with a radiolabelled compound (competition assays)</li> <li>Determining kinetics for association and dissociation of radiolabelled ligand to the receptor</li> <li>Determining protein and antibody affinity to membrane targets</li> <li>Evaluating allosteric binding</li> </ul>	<p>Receptor sources:</p> <ul style="list-style-type: none"> <li>Cell membranes (from wild-type or transfected cells)</li> <li>Tissue homogenates</li> <li>Whole cells</li> <li>Immobilized recombinant receptors</li> </ul> <p>Radioligands:</p> <ul style="list-style-type: none"> <li>[<sup>3</sup>H], [<sup>35</sup>S], [<sup>125</sup>I]-labelled compounds (commercially available or custom synthesized)</li> </ul>	<ul style="list-style-type: none"> <li><math>K_d</math></li> <li><math>B_{max}</math></li> <li><math>IC_{50}</math> and <math>K_i</math></li> <li><math>k_{on}</math></li> <li><math>k_{off}</math></li> </ul>	<ul style="list-style-type: none"> <li>High-sensitivity</li> <li>High-throughput</li> <li>Robust</li> <li>Works well on membrane receptors</li> <li>Tritium labelling can be used for small molecules without changing structure</li> <li><sup>125</sup>I custom labelling can be used for peptides, proteins, antibodies</li> </ul>	<ul style="list-style-type: none"> <li>Works best when radioligand <math>K_d &lt; \sim 50</math> nM</li> <li>A suitable high-affinity radioligand may need to be custom synthesized (if not commercially available)</li> <li>Cytoplasmic (soluble) proteins need to be first immobilized</li> </ul>
<b>Receptor Autoradiography</b>	Radiometric technique based on quantitative phosphor imaging of radioactivity in tissue sections from untreated or drug-treated animals (glass slides)	<ul style="list-style-type: none"> <li>Determining distribution of receptor-ligand binding sites</li> <li>Quantifying regional receptor levels</li> <li>Conducting <i>ex vivo</i> receptor occupancy studies to determine the percentage of receptor occupancy produced by a given dose of drug administered to the live animal</li> </ul>	<p>Receptor sources:</p> <ul style="list-style-type: none"> <li>Rodent tissue or organs from untreated or drug-treated animals</li> <li>Human tissues (post-mortem or from surgical removal)</li> <li>NHP tissues (post-mortem)</li> </ul> <p>Radioligands:</p> <ul style="list-style-type: none"> <li>[<sup>3</sup>H], [<sup>35</sup>S], [<sup>125</sup>I]-labelled compounds (commercially available or custom synthesized)</li> </ul>	<ul style="list-style-type: none"> <li><math>K_d</math></li> <li><math>B_{max}</math></li> <li>% Receptor occupancy versus dose</li> <li><math>ED_{50}</math></li> </ul>	<ul style="list-style-type: none"> <li>High-sensitivity</li> <li>Low amount of tissue required</li> <li>Can be used to detect binding to small tissue regions</li> <li>For <i>ex vivo</i> occupancy studies, lower cost than PET approaches and applicable to a wider range of receptors</li> <li>Tritium labelling can be used for small molecules without changing structure</li> <li><sup>125</sup>I custom labelling can be used for peptides, proteins, antibodies</li> </ul>	<ul style="list-style-type: none"> <li>Receptor should be present at sufficient abundance in the tissues</li> <li>A suitable high-affinity radioligand may need to be custom synthesized (if not commercially available)</li> <li><i>Ex vivo</i> studies may underestimate receptor occupancy if test (unlabelled) drug dissociates appreciably from the receptor during the <i>in vitro</i> incubation period in radioligand</li> </ul>
<b>Surface Plasmon Resonance (SPR)</b>	Label-free optical detection technique used to measure interacting biomolecules kinetics parameters (sensor chip)	<ul style="list-style-type: none"> <li>Testing antibody-antigen interaction</li> <li>Determining binding kinetics parameters of a variety of molecules (proteins, carbohydrates, small molecules etc.)</li> </ul>	<p>Ligands and analytes:</p> <ul style="list-style-type: none"> <li>Compounds</li> <li>Peptides</li> <li>Proteins and antibodies</li> <li>Saccharides</li> </ul> <p>Sample media:</p> <ul style="list-style-type: none"> <li>Buffer ± solvent</li> <li>Serum</li> <li>Plasma</li> </ul>	<ul style="list-style-type: none"> <li><math>K_d</math></li> <li><math>k_{on}</math></li> <li><math>k_{off}</math></li> <li>Selectivity/cross-reactivity</li> <li>Binding inhibition</li> <li>Analyte concentration</li> </ul>	<ul style="list-style-type: none"> <li>High-sensitivity</li> <li>High-reproducibility</li> <li>Label-free detection</li> <li>Low sample required</li> <li>Tolerance for sample impurity</li> <li>Real-time monitoring</li> </ul>	<ul style="list-style-type: none"> <li>Not suitable for whole cells</li> <li>Works best with non-membrane proteins</li> <li>Binding response is dependent on target and ligand MW ratio</li> <li>Small molecule binding to large receptor proteins hard to detect</li> </ul>
<b>ELISAs</b>	Immunoassay to detect target proteins using antibodies (96-well plates)	<ul style="list-style-type: none"> <li>Determining target protein expression levels</li> <li>Confirming cell surface expression of membrane receptors</li> <li>Assessing effect of drug treatments on the target expression profile</li> </ul>	<p>Target protein sources:</p> <ul style="list-style-type: none"> <li>Tissue homogenates</li> <li>Cell lysates</li> <li>Plasma</li> </ul>	<ul style="list-style-type: none"> <li>Target concentration</li> <li>Relative comparison of expression levels between tissues/treatments</li> </ul>	<ul style="list-style-type: none"> <li>High-sensitivity</li> <li>High-specificity</li> <li>Label-free detection</li> </ul>	<ul style="list-style-type: none"> <li>Lengthy protocols reduce high-throughput capabilities</li> <li>Relies on availability of antibodies that bind the target proteins</li> </ul>
<b>Fluorescence Assays</b>	Fluorescence-based assays used to measure downstream signalling of receptors upon ligand binding in live cells (96-well/384-well plates)	<ul style="list-style-type: none"> <li>Determining the effect of ligand-receptor interaction on cell signalling pathway(s)</li> <li>Profiling the pharmacology of the compound/receptors (agonist, inverse agonist, antagonist, allosteric modulators etc.)</li> </ul>	<ul style="list-style-type: none"> <li>Live cells (endogenous or transfected receptors)</li> </ul>	<ul style="list-style-type: none"> <li><math>E_{max}</math></li> <li><math>EC_{50}</math></li> <li><math>IC_{50}</math></li> </ul>	<ul style="list-style-type: none"> <li>High-throughput</li> <li>Suitable for functional studies</li> </ul>	<ul style="list-style-type: none"> <li>Some test molecules may exhibit intrinsic fluorescence (quenching) that might reduce accuracy and sensitivity</li> <li>Background fluorescence from cells may increase the signal-to-noise ratio</li> </ul>
<b>Cellular Uptake &amp; Release</b>	Radiometric technique used to measure transport of molecules into cells or release of labelled molecules from cells (96-well plate for uptake/perfusion chambers for release)	<ul style="list-style-type: none"> <li>Determining cellular uptake and release of labelled substrates, drugs, neurotransmitters or proteins</li> </ul>	<ul style="list-style-type: none"> <li>Live cells</li> <li>Synaptosomes</li> </ul> <p>Radioligands:</p> <ul style="list-style-type: none"> <li>[<sup>3</sup>H], [<sup>14</sup>C], [<sup>35</sup>S], [<sup>125</sup>I]-labelled compounds (commercially available or custom synthesized)</li> </ul>	<ul style="list-style-type: none"> <li>Compounds potency and efficacy against membrane transporters</li> <li>Compounds efflux from cells</li> <li><math>IC_{50}</math> and <math>K_i</math> for test compounds to inhibit cellular uptake of radiolabelled substrates</li> <li><math>V_{max}</math> and <math>K_m</math> for transport of labelled substrates</li> </ul>	<ul style="list-style-type: none"> <li>High-sensitivity</li> <li>Robust</li> <li>Uptake and release can be determined over time</li> <li>Tritium labelling can be used for small molecules without changing structure</li> <li><sup>125</sup>I custom labelling can be used for peptides, proteins, antibodies in uptake assays</li> </ul>	<ul style="list-style-type: none"> <li>Low throughput relative to radioligand binding assay approaches</li> <li>A suitable high-affinity radioligand may need to be custom synthesized (if not commercially available)</li> <li>Subcellular localization of labelled substrate cannot be determined without further processing</li> </ul>